medi GmbH & Co. KG
Medicusstraße 1
D-95448 Bayreuth
Germany
T +49 921 912-0
F +49 921 912-780
ortho@medi.de
www.medi.de

medi Australia Pty Ltd
83 Fennell Street
North Parramatta NSW 2151
Australia
T +61-2 9890 8696
F +61-2 9890 8439
sales@mediaustralia.com.au
www.mediaustralia.com.au

medi Austria GmbH
Adamgasse 16/7
6020 Innsbruck
Austria
T +43 512 57 95 15
F +43 512 57 95 15 45
vertrieb@medi-austria.at
www.medi-austria.at

medi Belgium NV
Staatsbaan 77/0099
3945 Ham
Belgium
T: +32-11 24 25 60
F: +32-11 24 25 64
info@medibelgium.be
www.medibelgium.be

medi Canada Inc / médi Canada Inc
597, Rue Duvernay, Verchères
QC Canada J0L 2R0
T: +1 450-583-3317 / +1 800 361 3153
F: +1 888-583-6827
service@medicanada.ca
www.medicanada.ca

MAXIS a.s.,
medi group company
Slezská 2127/13
120 00 Prague 2
Czech Republic
T: +420 571 633 510
F: +420 571 616 271
info@maxis-medica.com
www.maxis-medica.com

medi Danmark ApS
Vejlgaardsvvej 59
2665 Vallensbaek Strand
Denmark
T +45-70 25 56 10
F +45-70 25 56 20
kundeservice@sw.dk
www.medit danmark.dk

medi Bayreuth Espana SL
C/Canigo 2-6 bajos
Hospitalet de Llobregat
08901 Barcelona
Spain
T +34-932 60 04 00
F +34-932 60 23 14
medi@mediespana.com
www.mediespana.com

medi France
Z.I. Charles de Gaulle
25, rue Henri Farman
93297 Tremblay en France Cedex
France
T +33-1 48 61 76 10
F +33-1 49 63 33 05
infos@medi-france.com
www.medi-france.com

medi Hungary Kft.
Bokor u. 21.
1037 Budapest
Hungary
T +36 1 371-0090
F +36 1 371-0091
info@medi.hu
www.medi.hu

medi Japan K.K.
5-1 Nihonbashi Kabutocho Chuoku
103-0026 Tokyo
Japan
T: +81 3 6778 2590
F: +81 3 5847 7901
info@medi-japan.jp
www.medi-japan.co.jp

medi Nederland BV
Heusing 5
4817 ZB Breda
The Netherlands
T +31-76 57 22 555
F +31-76 57 22 565
info@medi.nl
www.medi.nl

medi

Manumed T

Handgelenkorthese mit Daumefixierung ·
Wrist brace with thumb support ·
Orthèse poignet pouce · Ortesis de
muñeca con soporte de pulgar

Gebrauchsanweisung. Instructions for use. Mode
d'emploi. Instrucciones de uso. Instruções para
aplicação. Istruzioni per l'uso. Gebruiksaanwijzing.
Brugsvejledning. Bruksanvisning. Návod k použití.
Uputa za uporabu. Инструкция по использованию.
Kullanma kilavuzu. Instrukcja załadania. Οδηγία
εφαρμογής. Інструкція з використання. Návod na
použitie. دليل الاستخدام. Instruțiuni de utilizare.
使用说明. 使用說明.

4 046938 573571

E007642 / 09.2020

medi. I feel better.

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Important notes

This medical device **MD** is made for single patient use only. If it is used for treating more than one patient, the manufacturer's product liability will become invalid. If undue pain or an unpleasant sensation develops while you are wearing the product, please consult your doctor or orthotist immediately. Do not wear the product over open wounds, and use it only as instructed by your doctor or orthotist.

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
medi Trading (Shanghai) Co. Ltd.
Room 1703, T1, CES WEST BUND
CENTER,
No. 277, Long Lan Road
Xuhui District
200232 Shanghai
People's Republic of China
T: +86-21 50582319
F: +86-21 50582319

Medi Turk Ortopedi Medikal İc ve Dis
Tic. AS
Mustafa Kemal Mah. 2156 Sokak No: 18
/4-5
06510 Çankaya Ankara
Turkey
T: +90 312 435 20 26
F: +90 312 434 17 67
info@medi-turk.com
www.medi-turk.com

medi UK Ltd.
Plough Lane
Hereford HR4 0EL
Great Britain
T +44-1432 37 35 00
F +44-1432 37 35 10
enquiries@mediuk.co.uk
www.mediuk.co.uk

medi Ukraine LLC
Evhena Sverstiuka str, 11.
Kiev 02002
Ukraine
T: +380 44 591 11 63
F: +380 44 392 73 73
info@medi.ua
http://medi.ua

medi USA L.P.
6481 Franz Warner Parkway
Whitsett, N.C. 27377-3000
USA
T +1-336 4 49 44 40
F +1-888 5 70 45 54
info@mediusa.com
www.mediusa.com

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Manumed T

Intended purpose

Manumed T is a wrist brace with thumb support used for immobilisation.

Indications

For all indications which require immobilisation of the wrist with thumb support, such as:

- Sprain on the wrist and/or thumb
- Post-operative immobilisation of the wrist and thumb
- Lesion of collateral ligament of the thumb basal joint
- Tendonitis, tendosynovitis
- Inflammatory phases for rheumatic conditions
- Trapeziometacarpal osteoarthritis

Contraindications

Unstable fractures in the area of the hand.

Risks / Side Effects

If the aids are too tight, this can cause pressure sores or constriction of blood vessels or nerves. Consult your doctor before using this device if you have the following conditions:

- A condition affecting the application area
- Paraesthesia, circulatory disorders
- Lymph drainage disorders

Tightly fitting aids can cause skin irritation (pressure, sweat, material).

Intended patient groups

Healthcare professionals should provide care to the adults and children, applying the available information on the measurements/sizes and necessary functions/indications and in line with the information provided by the manufacturer, and acting under their own responsibility.

Fitting instructions

- The premoulded aluminium splints can be taken out and easily moulded to different hand contours.
 - Open all belts and slip your hand into the brace (the aluminium splint is at the bottom - Fig. 1).
 - Close the front traction belt and then the back traction belt (Fig. 2).
 - Then close the broad middle traction belt. This belt is attached by hook and loop pads to the lower side of the brace and can be shifted before closing if necessary (Fig. 3).
 - Pull the thumb belt round the thumb (Fig. 4).
-

Care instructions

Please close the hook and loop fasteners before washing and remove the orthoses. Soap residues, lotions and ointments can cause skin irritation and material wear.

- Wash the product by hand, preferably using medi clean washing agent.
 - Do not bleach.
 - Leave to dry naturally.
 - Do not iron.
 - Do not dry clean.
-

Storage instructions

Keep the product in a dry place and do not expose to direct sunlight.

Material composition

Aluminum, polyamide.

Liability

The manufacturer's liability will become void if the product is not used as intended. Please also refer to the

corresponding safety information and instructions in this manual.

Disposal

The product can be disposed of in the domestic waste.

Your medi team

wishes you a quick return to full fitness.

In the event of any complaints regarding the product such as damage to the fabric or a fault in the fit, please report to your specialist medical retailer directly. Only serious incidents which could lead to a significant deterioration in health or to death are to be reported to the manufacturer or the relevant authorities in the EU member state. The criteria for serious incidents are defined in Section 2, No. 65 of the Regulation (EU) 2017/745 (MDR). The traceability of this product is assured via a UDI code UDI.

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