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## INSTRUCTIONS FOR USE GB

Dear Customer,

We would like to ask you to follow these Instructions for Use with care. If you should have any questions, please consult the doctor treating you or your nearest stockist, or contact us directly.

### FUNCTION

Specific fixation of the shoulder joint with free function and correction of the subluxation and medial rotation.

### INDICATIONS

- Hemiplegia
- Plexus injuries
- Skull/brain trauma
- Peripheral nerve damage

### CONTRA-INDICATIONS

Generally speaking you should see the doctor treating you about the indication and method of wearing an orthopaedic aid and if you have the following conditions:

- Allergic, inflammatory or lesion-specific skin alterations (e.g. swelling, reddening) of the body areas being treated.
  - Impaired circulation or lymphatic soft tissue swelling.
  - Neurogenically-specific organoleptic and skin trophic disorders in the body area being treated (sensory disorders with and without skin damage).
- So far no consequential incompatibilities or allergic reactions have been reported.

### SIDE-EFFECTS

To date, there have been no reports of serious general side effects with proper use and fitting.

Local pressure symptoms and impaired circulation can be prevented with sufficient certainty on an individual basis if allowance is made for any contraindications and with non-restrictive, consistently shaped shaped body fitting.

### FITTING

- 1. The shoulder cap A must be applied in such a manner that it fits comfortably on the shoulder.

**WARNING – Important information for the orthopaedic technician:** When the NEURO-LUX® Orthosis is applied for the first time it is necessary to ensure that the thermoplastic material of the shoulder cap A is carefully warmed with a hair dryer so that it can then be directly moulded onto the patient's shoulder.

- 2. The securing strap I is conducted on the opposite side underneath the armpit.

– 3. The arm cuff 2 must be applied in such a manner that the circular opening B lies over the elbow tip. The tensioning straps C and D are conducted in the direction of the shoulder joint.

- 4. The arm cuff then encircles the upper arm and forearm.

– 5. The closures must not be drawn too tight as this could result in blood stasis. This must be avoided!

- 6. Link the connecting straps C and D with the same coloured closures of the shoulder cap.

– 7. The straps must be lightly tensioned. The Orthosis must only be worn while awake because pressure sores can arise whilst asleep.

### IMPORTANT

The product is to be used by one patient only.

Improper modifications to the product and/or improper use of the product referred to above exempt the manufacturer from product liability.

You ought to discuss potential reciprocal health risks or other disadvantages with specific treatments that may arise in connection with the use of the product with the doctor treating you.

For the product to work effectively over a long period, it must not be worn in conjunction with fatty or acidic remedies, ointments or lotions.

### QUALITY MANAGEMENT SYSTEM

All SPORLASTIC GmbH products undergo product testing as part of our Quality Management System.

However, should you have any complaints about our product, please get in touch with your stockist or directly with us.

### CARE

We recommend that NEURO-LUX® be washed gently in lukewarm water with detergent for delicate items and that it then be air dried (do not dry on a heater).

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Wszystkie wyroby SPORLASTIC GmbH podlegają kontroli w ramach naszego systemu zarządzania jakością.