medi GmbH & Co. KG
Medicusstraße 1
95448 Bayreuth / Germany
T +49 921 912-0
F +49 921 912-780
ortho@medi.de
www.medi.de

medi Australia Pty Ltd
83 Fennell Street
North Parramatta NSW 2151
Australia
T +61 2 9890 8696
F +61 2 9890 8439
sales@mediaustralia.com.au
www.mediaustralia.com.au

medi Austria GmbH
Adamgasse 16/7
6020 Innsbruck
Austria
T +43-512 57 95 15
F +43-512 57 95 15 45
vertrieb@medi-austria.at
www.medi-austria.at

medi Belgium NV
Staatsbaan 77/0099
3945 Ham
Belgium
T: +32-11 24 25 60
F: +32-11 24 25 64
info@medibelgium.be
www.medibelgium.be

medi Canada Inc / médi Canada Inc
597, Rue Duvernay, Verchères
QC Canada J0L 2R0
T: +1 450-583-3317 / +1 800 361 3153
F: +1 888-583-6827
service@medicanada.ca
www.medicanada.ca

MAXIS a.s.,
medi group company
Slezská 2127/13
120 00 Prague 2
Czech Republic
T: +420 571 633 510
F: +420 571 616 271
info@maxis-medica.com
www.maxis-medica.com

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medi Norway AS
Vestheimvegen 35
4250 Kopervik
Norway
T: +47 52 84 45 00
F: +47 52 84 45 09
post@medinorway.no
www.medinorway.no

medi Polska Sp. z o.o.
ul. Łabędzka 22
44-121 Gliwice
Poland
T: +48-32 230 60 21
F: +48-32 202 87 56
info@medi-polska.pl
www.medi-polska.pl

medi Bayreuth Unipessoal Lda.
Rua do Centro Cultural, no. 43
1700-106 Lisbon
Portugal
T +351 21 843 71 60
F +351 21 847 08 33
medi.portugal@medibayreuth.pt
www.medi.pt

medi RUS LLC
Business Center NEO GEO
Butlerova Street 17
117342 Moscow
Russia
T +7-495 374 04 56
F +7-495 374 04 56
info@medirus.ru
www.medirus.ru

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Important notes

This medical device **MD** is made for single patient use only. If it is used for treating more than one patient, the manufacturer's product liability will become invalid. If undue pain or an unpleasant sensation develops while you are wearing the product, please consult your doctor or orthotist immediately. Do not wear the product over open wounds, and use it only as instructed by your doctor or orthotist.

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Stabimed®

Intended purpose

Stabimed is a knee support brace which limits flexion and extension.

Indications

All indications that require reinforcement of the physiological stabilisation of the knee joint and/or stress relief of the joint, such as:

- Following injury to the collateral ligament
- For knee instability
- Following meniscal injury
- For knee osteoarthritis

Contraindications

None known at present.

Risks / Side Effects

If the aids are too tight, this can cause local pressure sores or constriction of blood vessels or nerves. Therefore you should consult the doctor treating you before you use it if the following circumstances apply to you:

- Condition or injury of the skin in the area of application, particularly with infectious indications (excessive heat, swelling or reddening)
- Paraesthesia and circulatory disorders (e.g. as seen in diabetes, varicose veins)
- Lymph drainage disorders – equally, swelling of soft tissue adjacent to the area of application which is not related to a condition

When wearing snug-fitting aids, this can result in localised skin inflammation or irritation that is due to mechanical irritation of the skin (particularly in conjunction with sweating) or due to the composition of the material.

Intended patient groups

Healthcare professionals should provide care to the adults and children, applying

the available information on the measurements/sizes and necessary functions/indications and in line with the information provided by the manufacturer, and acting under their own responsibility.

Setting the flexion / extension limits Only to be carried out by an orthopedic technician.

Extension limit wedge, EXT (limits extension): 0°, 10°, 20°, 30°

Flexion limit wedge, FLEX (limits flexion): 0°, 10°, 20°, 30°, 45°, 60°, 75°, 90°

1. Only change the wedges if the orthoses are in the lateral supports.
2. 0° extension wedge is pre-installed in your knee brace.
3. To change the wedges, open the top and bottom latches. Now, simply remove the wedges.
4. In the wedge set, the replacement wedges are marked R – FLEX, L – FLEX or R – EXT, L – EXT and with the respective angles.
5. Insert the new flexion wedge first and the new extension wedge next into the joint. Then close the two latches to secure the wedges. Repeat the process with the other joint.
WARNING: The extension and flexion wedges always have to be inserted.
6. Move the brace joints to check if the wedges are securely seated.

Important notes

Only change the extension and flexion limits on instruction from your doctor or orthotist. In order to prevent hyperextension, always insert extension wedges or 0° wedges. The wedges define the exact excursion of the orthosis. In difficult soft-tissue conditions, e.g. where the soft-tissues are under compression, it may be necessary to insert a larger wedge to limit the joint movement to the desired excursion. In

addition, make sure that the orthosis joints are correctly positioned: the middles of the joints should be level with the upper border of the patella. The joint must lie behind the midline of the leg.

Fitting instructions

- Open all straps and the two hook and loop fasteners. Pull the two straps back a little so they don't hamper you when you put the brace on.
- Now put your leg into a 45-degree angle flexion position. Wrap the brace tightly around your leg from the back and close it with the two hook and loop strips (image 1).
- Make sure that the kneecap is in the middle of the patella recess.
- The centres of the brace joints should be at the top edge of the kneecap and just behind the lateral centre line of the leg (image 2).
- Now close the bottom strap and then the top strap (image 3).
- If possible, walk a few steps and check if the brace is properly fitted (image 4). Correct the fit of the brace at the leg, if required.

Care instructions

Please close the hook and loop fasteners before washing and remove the orthoses (image 5, 6, 7). Soap residues, lotions and ointments can cause skin irritation and material wear.

- Wash the product by hand, preferably using medi clean washing agent.
- Do not bleach.
- Leave to dry naturally.
- Do not iron.
- Do not dry clean.

Storage instructions

Keep the brace in a dry place and do not expose to direct sunlight.

Material composition

Avional, polyamide, elasthane, PU

Liability

The manufacturer's liability will become void if the product is not used as intended. Please also refer to the corresponding safety information and instructions in this manual.

Disposal

The product can be disposed of in the domestic waste.

Your medi team

wishes you a quick return to full fitness.

In the event of any complaints regarding the product such as damage to the fabric or a fault in the fit, please report to your specialist medical retailer directly. Only serious incidents which could lead to a significant deterioration in health or to death are to be reported to the manufacturer or the relevant authorities in the EU member state. The criteria for serious incidents are defined in Section 2, No. 65 of the Regulation (EU) 2017/745 (MDR). The traceability of this product is assured via a UDI code [UDI](#).

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